


Abstract

Liquid biopsy is a non-invasive sample source that can be utilized to assess cancer burden by measuring the tumor-derived fraction of circulating, cell-free DNA (cfDNA) from plasma. We evaluated two methods to monitor cancer burden using cfDNA: whole genome bisulfite sequencing (WGBS) and targeted amplicon sequencing for 56 oncology related genes. We also independently assessed whole genome sequencing of cell free DNA for use in other noninvasive assays. Genome-wide hypomethylation is a surrogate biomarker for cancer that can be detected independently of tumor genotype. Amplicon-based detection of tumor-specific mutations assesses tumorigenesis and potential therapeutic resistance. To compare the efficacy of both assays, cfDNA was extracted from tumor-bearing patients and normal controls. To monitor methylation density, WGBS was performed using 5 ng of bisulfite-converted cfDNA with the Accel-NGS[®] Methyl-Seq DNA Library Kit. To detect tumor-specific mutations, 10 ng of cfDNA was used for the Accel-Amplicon[™] 56G Oncology Panel. Six out of eight cancer samples demonstrated significant hypomethylation in cfDNA, ranging from 2 - 40% when compared to healthy controls. The 56 gene amplicon panel identified point mutations in the cfDNA of only three samples but which also had the highest observed hypomethylation (18-40%). For all but two cancer samples, corresponding mutations were also found in the primary tumor at allele frequencies significantly higher than in the cfDNA fraction (eg 22% in tumor vs. 5% in cfDNA). The three cancer samples that had primary tumor mutations that were not detected in cfDNA also had the lowest observed hypomethylation. Therefore, a correlation between hypomethylation and detection of tumor mutations in the cfDNA fraction may exist. For whole genome sequencing of cfDNA, the Accel-NGS 2S library preparation kit enabled PCR-free sequencing from inputs of 15 ng and produced uniform, comprehensive genome coverage, which provides a significant advantage for assessment of aneuploidy and copy number variation.


Superior Methyl-Seq Performance



1 ng Arabidopsis DNA Input	Swift Methyl-Seq	Traditional	Random Priming
% Reads Aligned	83.3%	80.7%	73.4%
Avg. Genome Coverage	18X	10X	12X
% Genome Covered ≥10x	77%	17%	31%
% Duplicate Reads	18%	62%	46%
Estimated Library Size	38 M	6 M	12 M

Paired-end sequencing was performed on a HiSeq with V4 chemistry with 125 bp PE. Analysis performed using BSMAP and Picard tools using 30M reads for each method for direct comparison.

Accel-Amplicon 56G Panel



56G Panel: Hotspot and Comprehensive Coverage of 56 Oncology-related Genes


Gene	Number of Amplicons	Hotspot Loci (White)	Comprehensive Coding Exon Coverage (Blue)
ABL1	5		
AKT1	2		
ALK	2		
APC	9		
ATM	19		
BRAF	2		
CDH1	3		
CDKN2A	2		
CSF1R	2		
CTNNA1	1		
DDR2	1		
DNMT3A	1		
EGFR	9		
ERBB2	4		
ERBB4	8		
EZH2	1		
FBXW7	6		
FGFR1	4		
FGFR2	4		
FGFR3	6		
FLT3	4		
FOXL2	1		
GNAI1	2		
GNAQ	2		
GNAS	2		
HNF1A	4		
HRAS	4		
IDH1	1		
IDH2	2		
JAK2	3		
JAK3	3		
KDR	9		
KIT	2		
KRAS	3		
MAP2K1	5		
MET	6		
MPL	1		
MSH6	4		
NOTCH1	3		
NPM1	14		
NRAS	3		
PDGFR4	5		
PIK3CA	6		
PTEN	14		
PTPN11	2		
RET	3		
SMAD4	3		
SMARCB1	4		
SMO	11		
SRC	14		
TP53	21		
TSC1	12		
VHL	6		
STK11	5		

56G panel genes and number of amplicons per gene. Hotspot loci (white), contiguous, overlapping coverage (blue) and comprehensive coding exon coverage for TP53 (darker blue).

Experimental Design

Tumor bearing blood, n = 8
(Streck Cell-free DNA BCT[®])

Healthy control blood, n = 5
(Streck Cell-free DNA BCT[®])




Bisulfite conversion and Accel-NGS Methyl-Seq library construction

Accel-Amplicon 56G library construction

Calculate hypomethylation status of cancer samples compared to healthy controls from 10M Illumina MiSeq[®] reads using MethyPipe.

LoFreq and GATK variant calling from 5000X coverage was compared to corresponding tumor FFPE and normal adjacent tissues.

Accel-NGS Methyl-Seq



- ✓ 2-hour library prep workflow
- ✓ Post-bisulfite prep uses efficient Adaptase ssDNA library technology
- ✓ Supports input range of 100 pg to 100 ng with minimal PCR cycles
- ✓ Retains high sequence complexity for comprehensive and uniform methylome coverage


Two Liquid Biopsy Assays: Genome-Wide Hypomethylation and Mutation Detection

mL Plasma	ng/mL cfDNA	Pathology	cfDNA Hypomethylation	56G Mutation	Normal Adjacent	FFPE Tumor	cfDNA
2.5	6.3	Fallopian tube high-grade papillary serous carcinoma	0.4 %	TP53 E285K	0%	48%	0%
5.0	4.3	5 cm ovarian 'borderline' serous content	1.1 %	BRAF V600E	0%	14%	0%
3.8	4.4	Recurrent pT2, pN0 mammary carcinoma	2.4 %	PIK3CA H1047R	0%	17%	0%
4.0	10.5	pT1/pN1 pancreatic adenocarcinoma	3.6 %	-	-	-	-
3.0	6.7	Metastatic colon cancer to the liver	4.4 %	-	-	-	-
4.5	7.1	14 cm ovarian 'borderline' serous content	18.0 %	BRAF V600E	0%	23%	1%
4.5	2.6	Colon-cancer, non-resectable Adenocarcinoma	18.0 %	TP53 frameshift exon 8	0%	15%	2%
4.5	2.9	Metastatic colorectal adenocarcinoma	43.4 %	PIK3CA E545K APC Q1429* TP53 Q38* KRAS G13D	0% 0% 0% 0%	23% 20% 21% 22%	11% 5% 14% 5%

A correlation between cfDNA hypomethylation and detection of tumor mutations in cfDNA may exist. Significant hypomethylation was detected in 6 of the 8 samples, and the 56 gene amplicon panel identified point mutations in the cfDNA of the three samples with the highest observed hypomethylation. Concordance was observed between corresponding cfDNA and FFPE tumors, when mutations were detected.

Limit of Detection (gDNA)


1 ng of gDNA = 334 chromosomal copies of any locus
Achieving 1% detection of a mutation → 3 chromosomal copies



To assess the limit of detection of Accel-NGS 2S Hyb, DNA samples from two individuals (A and B) with different ethnic backgrounds were used to prepare libraries. 100 ng of DNA from one individual with a 0.5% or 1% spike-in of the DNA from the second individual was used as the input DNA. Once libraries were prepared, they were hybridized to xGen Pan-Cancer probes and SNPs were detected within this panel.


chr: POS	ALLELE FREQUENCIES				
	A Back-ground	B Back-ground	C Back-ground	1% A into 10 ng B	1% A into 10 ng C
2: 212244718	100%	0%	0%	0.6%	1.0%
12: 25361074	100%	0%	0%	1.6%	1.9%
12: 25361142	100%	0%	0%	1.1%	0.9%
12: 25361646	100%	0%	0%	1.9%	1.6%
12: 40688695	100%	0%	0%	0.5%	1.1%
12: 115108136	100%	0%	0%	0.7%	2.0%

To determine if SNPs present at 1% allele frequency could be detected, 1% of cfDNA sample (A) with a unique ethnic background was spiked into two, 10 ng cfDNA samples (B and C) of different ethnic backgrounds. SNPs at 100% from sample A could be detected around 1% when those SNPs are not present in B and C backgrounds. Libraries were sequenced to an average coverage of 8700X.



The circos plot depicts methylation density of 1 Mb bins across chromosomes 1-22 for the metastatic colorectal adenocarcinoma sample, where red is hypomethylated (>3 SD lower than normal mean MD) and green is comparable to normal.

Accel-NGS 2S Hyb



- ✓ Simple with-bead protocol
- ✓ Broad input range 10 pg - 1 µg
- ✓ Repairs both 5' and 3' termini to enhance ligation efficiency
- ✓ Compatible with cfDNA and FFPE samples
- ✓ Increased library complexity
- ✓ Balanced coverage of AT-/ GC-rich regions

Conclusions

Accel-NGS Methyl-Seq library prep kit:
- Provides uniform, comprehensive methylome coverage.
- Enables liquid biopsy for genome-wide hypomethylation from 5 ng cfDNA.

Accel-Amplicon 56G targeted sequencing panel:
- Provides quality performance with > 95% on target and > 95% coverage uniformity.
- Enables a limit of mutation detection of 1% for liquid biopsy from 10 ng cfDNA.

Using both kits for liquid biopsy, a correlation was observed between percent hypomethylation and mutation detection for the tumor bearing cfDNA sample set presented.

The Accel-NGS 2S library prep kit:
- Provides uniform, comprehensive genome coverage.
- Enables PCR-free cfDNA sequencing from 10-15 ng input.

Thanks to Brian Gerwe, Alex Lopez and Mike Benway (Perkin Elmer) for cfDNA extraction and sequencing