



## ***Run Setup and Bioinformatic Analysis***

**Accel-NGS<sup>®</sup> 2S MID Indexing Kits**

# Sequencing MID Libraries

For MiSeq<sup>®</sup>, HiSeq<sup>®</sup>, and NextSeq<sup>®</sup> instruments:

- Modify the config file to create a fastq for index reads
- Using the Illumina<sup>®</sup> Experiment Manager software, specify 2 index reads for the run.
- In the CSV file, specify the MID index and the sample index sequences.
  - Samples will be demultiplexed based only on their sample index.
- Using a custom script, join MIDs to their respective fastq read headers, align these fastq, and analyze the reads with a common genomic coordinate.
- **The following slides contain instructions for the MiSeq instrument, but these can also apply to HiSeq and NextSeq instruments as well.**

# The Swift MID Adapter Sequence

Below is the sequence of the standard LT P7 adapter:

GATCGGAAGAGCACACGTCTGAACTCCAGTCAC[i7]ATCTCGTATGCCGTCTTCTGCTTG

Below is the sequence of the Swift MID P5 adapter, where 9 random N bases (instead of the standard Index 2, 8 bp D501-508 index sequence) is inserted:

AATGATACGGCGACCACCGAGATCTACACNNNNNNNNNACACTCTTTCCCTACACGACGCTCTTCCGATCT

# Pre-Run Instrument Set-Up

## Modify the config file to allow generation of an index fastq file during data analysis:

1. Stop the MiSeq<sup>®</sup> Reporter process (if it is running).
2. Locate the “MiSeq Reporter.exe.config” file located in C:/Illumina/MiSeq Reporter.
3. Open config file and search for a line that reads:  
“<add key=“CreateFastqForIndexReads” value=“0”/>”.
  - If this line is present, change the value from “0” to “1”.
  - If this line is not present, add the line to the config file using the add keys function under the app Settings tab.
4. Restart the MiSeq reporter process.
5. Re-queue the run for data analysis (if required).

# Instrument Set-Up and Sample Sheet Preparation

- Experimental set-up on Illumina<sup>®</sup> Experiment Manager software
  - Select MiSeq
  - On the MiSeq Application Selection page, select Other --> FASTQ Only
  - Enter Reagent Cartridge Barcode, choose TruSeq HT, choose “2” for Index Reads, ... Read Type: Paired End, enter desired number of cycles for each read, and uncheck all the boxes on the right (including adapter trimming)
    - **By selecting HT, we can alter the index sequences read in the sample sheet**
    - **2 Index reads ensures P7 index (Index 1) is read, as well as the P5 MID located in place of the P5 index (Index 2)**

The screenshot displays two panels from the Illumina Experiment Manager software. The left panel, titled "FASTQ Only Run Settings", contains the following fields and options:

- Reagent Cartridge Barcode\*: [Empty text box]
- Sample Prep Kit: [TruSeq HT dropdown menu]
- Index Reads: [Radio buttons for 0, 1, and 2; 2 is selected]
- Experiment Name: [Empty text box]
- Investigator Name: [Empty text box]
- Description: [Empty text box]
- Date: [12/ 4/2015 date picker]
- Read Type: [Radio buttons for Paired End and Single Read; Paired End is selected]
- Cycles Read 1: [151 spinner box]
- Cycles Read 2: [151 spinner box]

The right panel, titled "FASTQ Only Workflow-Specific Settings", contains the following checkboxes, all of which are unchecked:

- Custom Primer for Read 1
- Custom Primer for Index
- Custom Primer for Read 2
- Reverse Complement
- Use Adapter Trimming
- Use Adapter Trimming Read 2

\* - required field

# Instrument Set-Up and Sample Sheet Preparation

- On the next step, enter your index sequence if using high throughput indices
- If using low throughput indices, select a random index sequence until sample sheet status is valid
  - We will specify the low throughput sequence in the sample sheet
- For I5 Sequence, enter a random index number on the pull down menu
  - We will specify the MID (NNNNNNNNNN) in the sample sheet
  - Click Finish to generate the CSV file

Sample Sheet Wizard - Sample Selection

Samples to include in sample sheet \*-required field  Maximize

Sample ID*	Sample Name	Plate	Well	Index1 (I7)*	I7 Sequence	Index2 (I5)*	I5 Sequence	Sample Project	Description
1				D701	ATTACTCG	D501	TATAGCCT		
2				D702	TCCGGAGA	D501	TATAGCCT		
3				D703	CGCTCATT	D503	CCTATCCT		
4				D704	GAGATTCC	D504	GGCTCTGA		

?

Sample Sheet Status: Valid  
Reason:

# Instrument Set-Up and Sample Sheet Preparation

- Alter your sample sheet to represent the real index sequences
  - When using a low throughput index, include the index and the next 2 bases
    - » For example with I2: Enter **CGATGTAT**. Bold represents the index whereas the underlined bases represent the next two bases on the adapter
  - For Index 2, enter the mID sequence: NNNNNNNNN
  - During the MiSeq run, the samples will be separated based only on their Index 1 indexes, since any Index 2 reads will be valid

	A	B	C	D	E	F	G	H	I	J
1	[Header]									
2	IEMFileVersion		4							
3	Investigator Name	Dr. Lenhart								
4	Experiment Name	mID								
5	Date	12/2/2015								
6	Workflow	GenerateFASTQ								
7	Application	FASTQ Only								
8	Assay	TruSeq HT								
9	Description									
10	Chemistry	Amplicon								
11										
12	[Reads]									
13		76								
14		76								
15										
16	[Settings]									
17	ReverseComplement		0							
18										
19	[Data]									
20	Sample_ID	Sample_Name	Sample_Plate	Sample_Well	I7_Index_ID	index	I5_Index_ID	index2	Sample_Project	Description
21		1 Sample #1			D701	CGATGTAT	D501	NNNNNNNNN		
22		2 Sample #2			D702	CAGATCAT	D502	NNNNNNNNN		
23		3 Sample #3			D703	TGACCAAT	D503	NNNNNNNNN		
24		4 Sample #4			D704	GCCAATAT	D504	NNNNNNNNN		

# Instrument Set-Up and Sample Sheet Preparation

- After the run, all reads will be separated based on Index 1 reads

Analysis	Imaging	Summary	Tile Status	TruSeq Controls	Indexing
<b>Reads mapped to Index Id</b>					
<b>Total Reads</b>	<b>PF Reads</b>	<b>% Reads Identified (PF)</b>			
30804952	27330618	98.0736			
<b>Index Number</b>	<b>Sample Id</b>	<b>Project</b>	<b>Index 1 (I7)</b>	<b>Index 2 (I5)</b>	
1	1	NA	CGATGTAT	NNNNNNNNN	
2	2	NA	CAGATCAT	NNNNNNNNN	
3	3	NA	TGACCAAT	NNNNNNNNN	
4	4	NA	GCCAATAT	NNNNNNNNN	

Non-Identified (Undetermined) reads are due to either poor quality or reads containing absent index sequences.



# Options for Retrieving Index Sequence Files

Instrument	bcl2fastq	MiSeq reporter	BaseSpace
MiniSeq	✓		
MiSeq	✓	✓	
NextSeq 500	✓		
HiSeq 2500	✓		
HiSeq 4000	✓		

If for any reason the data is not extracted correctly by the Sample Sheet setup, [bcl2fastq Conversion Software](#) can be used to correctly extract the data.

# Logic for Bioinformatic Analysis of ChIP-Seq and cfDNA

- After the sequencing run, all reads will be separated based on Index 1 only. Non-identified (or undetermined) reads are either due to poor quality or reads missing index sequences.
- Using a custom script, join the MIDs to the respective fastq read headers.
- Align these fastq using BWA.
- Using another custom script, analyze the reads with a common genomic coordinate.
  - If the reads have unique MID sequences, they represent fragment duplicates and should be retained as unique reads
  - If the reads have identical MID sequences, they represent PCR duplicates. Mark all reads except one as duplicates based on mapping and base quality.

# Logic for Bioinformatic Analysis of Low Frequency Variants

- After the sequencing run, all reads will be separated based on Index 1 only. Non-identified (or undetermined) reads are either due to poor quality or reads missing index sequences.
- Using a custom script, join the MIDs to the respective Fastq read headers.
- Align these Fastq using BWA.
- Use Picard Tools or Samtools to collect metrics like % duplication, reads on target, etc.
- Using another custom script, group fragments with at least 3 PCR duplicate reads per MID for variant calling (also considering Flag and Cigar values). Within the group, determine a consensus sequence for each fragment which eliminates sequencing and PCR errors present at less than 50%.
  - The following publication contains some details concerning this kind of analysis, and may be a useful reference:
    - S.R. Kennedy, et al. Nat Protoc. 2014 November; 9(11):2586-2606

# Fulcrum Genomics Analysis Tools

- To perform this analysis, we recommend using UMI tools from the fgbio package (<https://github.com/fulcrumgenomics/fgbio>)
  - It enables you to tag your bam file with the MIDs from your Index 2 fastq (AnnotateBamWithUmis)
  - Use the Picard MarkDuplicate with “BARCODE\_TAG=RX” option to get duplicate rate based on MIDs.
  - Use “GroupReadsByUmi” to group reads with the same MID.
  - Use “CallMolecularConsensusReads” to enrich for MID reads for variant calling.



# THANK YOU

[www.swiftbiosci.com](http://www.swiftbiosci.com)

© 2016, Swift Biosciences, Inc. The Swift logo is a trademark and Accel-NGS is a registered trademark of Swift Biosciences. Illumina, TruSeq, MiSeq, HiSeq, and NextSeq are registered trademarks of Illumina, Inc. 16-1188, 10/16

**Swift**   
BIOSCIENCES™