

Set S3 Adapters	Sequence	26796
Reagent Y2 (Index S749)	TCCAGTCG	24 µl
Reagent Y2 (Index S750)	TGTATGCG	24 µl
Reagent Y2 (Index S751)	TCATTGAG	24 µl
Reagent Y2 (Index S752)	TGGCTCAG	24 µl
Reagent Y2 (Index S753)	TATGCCAG	24 µl
Reagent Y2 (Index S754)	TCAGATTC	24 µl
Reagent Y2 (Index S755)	GGTTGGAC	24 µl
Reagent Y2 (Index S756)	GACACTTA	24 µl
Reagent Y2 (Index S757)	GCTATGGA	24 µl
Reagent Y2 (Index S758)	GTAACCGA	24 µl
Reagent Y2 (Index S759)	GGCAAGCA	24 µl
Reagent Y2 (Index S760)	GAACGACA	24 µl
Reagent Y2 (Index S761)	GCGTCGAA	24 µl
Reagent Y2 (Index S762)	AAGGCGAT	24 µl
Reagent Y2 (Index S763)	CAGGCATT	24 µl
Reagent Y2 (Index S764)	AACTGTAT	24 µl
Reagent Y2 (Index S765)	ATGCTTGA	24 µl
Reagent Y2 (Index S766)	AGTATCTG	24 µl
Reagent Y2 (Index S767)	ATGTAATG	24 µl
Reagent Y2 (Index S768)	ACACATGT	24 µl
Reagent Y2 (Index S769)	ATAGCACG	24 µl
Reagent Y2 (Index S770)	ATATTGTA	24 µl
Reagent Y2 (Index S771)	CAATTGAT	24 µl
Reagent Y2 (Index S772)	CACGTCGT	24 µl

Set S4 Adapters	Sequence	26896
Reagent Y2 (Index S773)	AGTCTGTA	24 µl
Reagent Y2 (Index S774)	CCGTATCT	24 µl
Reagent Y2 (Index S775)	CGTCTCCT	24 µl
Reagent Y2 (Index S776)	CAAGACCT	24 µl
Reagent Y2 (Index S777)	CCTAGTAT	24 µl
Reagent Y2 (Index S778)	CCACCGAT	24 µl
Reagent Y2 (Index S779)	CTATCATG	24 µl
Reagent Y2 (Index S780)	CATGAATG	24 µl
Reagent Y2 (Index S781)	CTGTACGG	24 µl
Reagent Y2 (Index S782)	CACTCGAG	24 µl

Set S4 Adapters	Sequence	26896
Reagent Y2 (Index S783)	CCGACAAG	24 µl
Reagent Y2 (Index S784)	CTTGCTTC	24 µl
Reagent Y2 (Index S785)	CGCCTTAT	24 µl
Reagent Y2 (Index S786)	GCAACCAT	24 µl
Reagent Y2 (Index S787)	TGACCGTT	24 µl
Reagent Y2 (Index S788)	TTGAGCTC	24 µl
Reagent Y2 (Index S789)	CCACATTG	24 µl
Reagent Y2 (Index S790)	AGCCAACCT	24 µl
Reagent Y2 (Index S791)	ATCACGTT	24 µl
Reagent Y2 (Index S792)	TCTCGGTT	24 µl
Reagent Y2 (Index S793)	TTGACTCT	24 µl
Reagent Y2 (Index S794)	TCGAAGTG	24 µl
Reagent Y2 (Index S795)	CACCCAAA	24 µl
Reagent Y2 (Index S796)	CTTCACAT	24 µl

During library prep, make sure to note which indexed adapter you are using with your sample and do not use the same indexed adapter on two different samples you plan to multiplex together.

Reagents	26596	26696	26796	26896	269384
Reagent B2	211 µl	211 µl	211 µl	211 µl	844 µl
Reagent R1	528 µl	528 µl	528 µl	528 µl	2,112 µl

For questions, please contact TechSupport@swiftbiosci.com.

Section C: Helpful Information and Troubleshooting

Problem	Possible Cause	Suggested Remedy
Library migrates unexpectedly on Bioanalyzer.	When analyzed on the Agilent High Sensitivity chip, migration behavior overestimates library size of PCR-free libraries made from DNA fragmented to the 200–300 base range (as required in this protocol).	<ul style="list-style-type: none">• “200 bp insert” PCR-free libraries should migrate to a ≈500 bp peak and “350 bp insert” PCR-free libraries should migrate to a ≈800 bp peak on the High Sensitivity Chip.• Consult the Expected Results section and the application note released by Covaris titled “Analysis of DNA Fragments Using the Agilent 2100 Bioanalyzer” to ensure proper analysis of library size.
	Over-amplification of library leads to the formation of heteroduplex structures that migrate abnormally.	<ul style="list-style-type: none">• Quantify library by qPCR, as other quantification methods will not accurately detect heteroduplex library molecules.• Perform the minimum number of PCR cycles necessary to avoid over-amplification.
DNA does not fragment properly: broad or top-sided (high molecular weight) sonication profile of fragmented DNA.	Impure DNA or fragmentation device malfunction.	<ul style="list-style-type: none">• Isopropanol purification, bead clean-up, column purification, or other method before fragmentation.• Ensure fragmentation device is functioning within manufacturer’s parameters.
Incomplete resuspension of beads after ethanol wash during SPRI™ steps.	Over-drying of beads.	Continue pipetting the liquid over the beads to break up clumps for complete resuspension.
Shortage of enzyme reagents.	Pipetting enzymes at -20 °C instead of 0-4 °C.	Allow enzyme reagents to equilibrate to 0-4 °C for 10 minutes prior to pipetting.
Retention of liquid in pipette tip	Viscous reagents may stick to pipette tip, especially for non-low retention tips.	Pipette up and down several times to ensure all liquid and/or beads are released from the pipette tip.

If you experience problems with your library prep, please contact us at TechSupport@swiftbiosci.com, or by phone at 734.330.2568 (9:00 am-5:00 pm ET, Monday-Friday).

General Warranty

Swift Biosciences, Inc. ("Swift") warrants that its products meet Swift's specifications at the time of delivery. Any sample or model used in connection with Swift's product literature is for illustrative purposes only and does not constitute a warranty that the products will conform to the sample or model.

To the maximum extent permitted by applicable law, Swift hereby expressly disclaims, and the buyer hereby expressly waives, any warranty regarding results obtained through the use of the products including, without limitation, any claim of inaccurate, invalid, or incomplete results. All other warranties, representations, terms and conditions (statutory, express, implied or otherwise) as to quality, condition, description, merchantability, fitness for purpose, or non-infringement (except for the implied warranty of title) are hereby expressly excluded.

All warranty claims on products must be made in writing within ninety (90) days of receipt of the products. Swift's sole liability and the buyer's exclusive remedy for a breach of this warranty is limited to replacement or refund at the sole option of Swift.

The warranties identified in this paragraph are Swift's sole and exclusive warranties with respect to the products and are in lieu of all other warranties, statutory, express or implied, all of which other warranties are expressly disclaimed, including without limitation any implied warranty of merchantability, fitness for a particular purpose, non-infringement, or regarding results obtained through the use of any product (including, without limitation, any claim of inaccurate, invalid or incomplete results), whether arising from a statute or otherwise in law or from a course of performance, dealing or usage of trade.

Limitation of Liability

Swift Biosciences, Inc. ("Swift") shall have no liability under the warranties cited above with respect to any defect in the products arising from: (i) specifications or materials supplied by the buyer; (ii) wilful damage or negligence of the buyer or its employees or agents; (iii) abnormal working conditions at the buyer's premises; (iv) failure to follow Swift's use restrictions or instructions (whether oral or in writing); (v) misuse or alteration of the products without Swift's approval; or (vi) if the buyer is in breach of its payment obligations in regards to purchasing the products.

To the fullest extent allowed by law, in no event shall Swift be liable, whether in contract, tort, strict liability, negligence, warranty, or under any statute or on any other basis for any special, incidental, indirect, exemplary, punitive, multiple or consequential damages sustained by the buyer or any other person or entity arising out of or caused by product, Swift's performance or failure to perform its obligations relating to the purchase of product or performance of services, Swift's breach of these terms, the possession or use of any product, or the performance by Swift of any services, whether or not foreseeable and whether or not Swift is advised of the possibility of such damages, including without limitation damages arising from or related to loss of use, loss of data, downtime, procurement of substitute products or services, or for loss of revenue, profits, goodwill, or business or other financial loss.

The total liability of Swift arising under or in connection with the purchase of the products, including for any breach of contractual obligations and/or any misrepresentation, misstatement or tortious act or omission (including without limitation, negligence and liability for infringement of any third party intellectual property rights) shall be limited to damages in an amount equal to the amount paid to Swift under the purchase agreement.

The exclusion of liability shall apply only to the extent not prohibited by applicable law.

Notice to Purchaser: Limited License

This product is for research use only and is licensed to the user under Swift Biosciences intellectual property only for the purchaser's internal purposes. Not for use in diagnostic procedures.



Swift Biosciences, Inc.

674 S. Wagner Road, Suite 100 • Ann Arbor, MI 48103 • 734.330.2568 • www.swiftbiosci.com

©2017, Swift Biosciences, Inc. The Swift logo is a trademark and Accel-NGS is a registered trademark of Swift Biosciences. This product is for Research Use Only. Not for use in diagnostic procedures. Illumina, HiSeq, MiniSeq, MiSeq, NextSeq, and TruSeq are registered trademarks of Illumina, Inc. Agencourt and AM-Pure are registered trademarks and SPRI, SPRIPlate, and SPRIselect are trademarks of Beckman Coulter, Inc. Qubit is a registered trademark and DynaMag is a trademark of Thermo Fisher Scientific, Inc. GenomePlex is a registered trademark and SeqPlex is a trademark of Sigma-Aldrich. Covaris is a registered trademark of Covaris, Inc. ITaq is a trademark of Bio-Rad Laboratories, Inc. Oligonucleotide sequences
© 2016 Illumina, Inc. All rights reserved. 17-1640, 08/17